


Design Collection



| | |
|----------------------------------|--------|
| Profondità seduta | 520 mm |
| Larghezza seduta | 530 mm |
| Altezza schienale (dalla seduta) | 68 cm |
| Lunghezza totale aperta | 168 cm |
| Altezza poggiatesta | 60 cm |
| Altezza schienale da aperta | 85 cm |
| Portata massima | 150 kg |

Colori



Verde Marrone

2 Lift

Poltrona GILDA ELEVATION è costituita da una struttura in legno e da meccanica in acciaio, rivestita in microfibra effetto mano-pesca anallergica. È dotata di 4 ruote piroettanti, di cui due con freno di stazionamento, maniglione di spinta e cintura di contenimento.

Facile estrazione dei braccioli e delle alette poggiatesta per consentire un facile passaggio dalla poltrona al letto o alla sedia a rotelle. La poltrona è dotata di 2 motori (sistema 1+1): un motore consente i classici movimenti letto e lift-up, mentre il secondo è per il sollevamento verticale con un'elevazione della poltrona fino ad un massimo di 16 cm.

The GILDA ELEVATION armchair is made of wood and mechanized steel structure, coated in micro-fiber hand-fishing anallergic effect. It is equipped with 4 swivel castors, two of which have a parking brake, a push handle and a seat belt.

Easy extraction of the armrests and headrest wings to allow easy passage from the armchair to the bed or to the wheelchair. The armchair is equipped with 2 motors (1 + 1 system): an engine allows classic lift and lift movements, while the second is for vertical lifting with an armchair up to a maximum of 16 cm.

Posizione relax



Posizione lift-up



Posizione elevation

